







Wet Lab Look-in Presentation

Ecophysiology and “Biosensing”: What can invertebrates tell us about water quality?

Sam Gurr



Lighthouse Point
New Haven, CT



Ecophysiology



Stony Brook
University




 AUTODESK®

The stressor: Hypoxia

Hypoxia

- Low dissolved oxygen (DO)

demand > supply



The stressor: Hypoxia


Hypoxia

- Periods of low dissolved oxygen (DO)

demand > supply

Diel-cycling hypoxia

- Varies on a seasonal and daily time scale



Common effects of diel-cycling hypoxia

- mass mortality
- growth, calcification
- behavior
- early life stage development
- calcification
- immunoregulatory response

LIMITATION!


Integrate over relatively
long time scales

Physiological “windows” of optimal aerobic performance

Aerobic performance

- Respiration rate
- Cardiac activity

This window can be **narrowed** by environmental stressors



Respiration rate of marine invertebrates:

Responses to oxygen decline and hypoxia

Conformer
(oxyconformity)


VS

Regulator
(oxyregulatory)

- Linear dependence on ambient oxygen availability

- Non-linear response
- Hits a clear transitional breakpoint at P_c

Critical P_{O_2} or P_c
onset of anaerobic metabolism



Rational

To understand and identify effects of *dynamic environmental stressors*...


Need a dynamic response!!!

Heartbeat rate



Non-invasive
Infrared
sensors

Important **sub-lethal**
physiological implications
for whole animal
metabolism



Rational

To understand and identify effects of *dynamic environmental stressors*...

Need a dynamic response!!!

Heartbeat rate



Non-invasive
Infrared
sensors

Common effects of diel-cycling hypoxia

- mass mortality
- growth, calcification
- behavior
- early life stage development
- calcification
- immunoregulatory response



Important **sub-lethal**
physiological implications
for whole animal
metabolism

Methods



1. *A. irradians* alters cardiac activity under exposure to *in-situ* diel-cycling dissolved oxygen


- 8 *A. irradians* per site
- 60 seconds of heartbeat data for each individual, cycle **repeated every 10 minutes**
- DO and temperature recorded every 15 minutes, Long Island Water Quality Index program




In-situ Heartbeat Deployments

- *A. irradians* heartbeat (bpm)
- Dissolved oxygen (mg L^{-1})
- Temperature ($^{\circ}\text{C}$)


Fire Island





Nicoll Bay



Seatuck



In-situ Cardiac Activity: Normoxic VS. Hypoxic conditions





In-situ Cardiac Activity: Normoxic VS. Hypoxic conditions

➡ Photosynthesis

➡ Organic matter mineralization


- **Heartbeat rate:**
dynamic daily
range of 20 bpm



In-situ Cardiac Activity: **Diel-cycling hypoxia**


- Daily oscillations of activity are **not temperature driven**

• *A. irradians* heartbeat (bpm)
• Dissolved oxygen (mg L^{-1})
• Temperature ($^{\circ}\text{C}$)




In-situ Cardiac Activity: **Diel-cycling hypoxia**

- Cardiac oscillations are **driven by DO decline**





In-situ Cardiac Activity: **Diel-cycling hypoxia**



- Cardiac activity always peaked when DO decline to **5 mg L⁻¹** during early to late mornings
- Evidence of a potential onset of:
 - decline of aerobic function
 - transition to anaerobic respiration



Cardiac response of diel-cycling dissolved oxygen (CRD_{DO})



Cardiac response of diel-cycling dissolved oxygen (CRD_{DO})



Cardiac response of diel-cycling dissolved oxygen (CRD_{DO})

Phase 1



Heartbeat
rate




Dissolved
oxygen

Oxyregulatory response

- maintain oxygen uptake and aerobic metabolism as DO becomes less available

Heartbeat rate change: +10 bpm

Duration: 8 – 10 hours (longest phase)



Cardiac response of diel-cycling dissolved oxygen (CRD_{DO})

Phase 2



Heartbeat
rate



Dissolved
oxygen


Transition to oxyconformity

- **Peak heartbeat rate at 5 mg L⁻¹**
- May indicate an initiation of anaerobic pathways

Remember: Pc, TcII, and ABT

Heartbeat rate change: -10 bpm

Duration: 4 – 4.5 hours (shortest phase)



Cardiac response of diel-cycling dissolved oxygen (CRD_{DO})

Phase 3



Heartbeat
rate




Dissolved
oxygen

"Stress and rest"

- Cardiac activity continues to decline to a minimum rate although DO increases
- **Minimum heartbeat rate at 5 – 7 mg L⁻¹**

Heartbeat rate change: -10 bpm

Duration: 5 – 6 hours



Cardiac response of diel-cycling dissolved oxygen (CRD_{DO})

Phase 4



Heartbeat
rate




Dissolved
oxygen

Cardiac and aerobic recovery



- Only phase when both heartbeat rate and DO **increase**
- Suggests an initial effort to restore aerobic function to basal heartbeat rates

Heartbeat rate change: +10 bpm

Duration: 4 – 6 hours




In a metaphorical sense...



As if these scallops reached the summit of Mt. Everest
every morning






THE RIVER PROJECT

Project Goals

- Coastal water quality determined by...



Bio-sensors

organisms used to detect complex conditions

- Goal:** use biosensors as an augmentation of traditional digital sensor data to better understand ecosystem health


How? Cardiac activity and valve gape behavior

- Blue mussel**
 - Atlantic coast native
 - filters water to feed and respire
 - semi-sessile (sedentary)






Process and Methods

Site collection


-  = New York City
-  = Stony Brook



- **Valve gape**
Hall effect sensor
- **Heartbeat rate**
Infrared sensor




Results



- Complex coupled dynamics from digital sensors
- Obvious influence of semidiurnal tides
- Interplay between
fresh water (Hudson River)
+
salt water (Atlantic Ocean)


Results

■ Digital sensors




Time

■ Biosensors



Results




Local

Tranlocated
(naïve)

Results

What does this mean?



In Summary

- Bio-sensors give a **unique view** of ecosystem status unachievable by digital sensors alone

Questions?

GAPE! QUICK!!!
That way he can't apply
a heartbeat sensor and
deploy us to our doom!

ERRRGH! Yup.
...he glued me shut.

